


High-Resolution Spectroscopy for AMO and Materials Science • Beamline 9.3.1

Berkeley Lab • University of California

Beamline Specifications

Photon Energy Range (keV)	Photon Flux (photons/s)	Spectral Resolution (E/ΔE)	Spot Size (mm)	Availability
2.2–6.0 [with Si(111) crystals]	>10 ¹¹	3000–8000	<0.5	NOW




Schematic layout of Beamline 9.3.1.

Beamline 9.3.1 is a bend-magnet beamline that provides high-resolution x-ray beams for electron, ion, and x-ray spectroscopy measurements in atomic, molecular, and materials science. The beamline is built around a double-crystal monochromator to provide higher-energy x-rays than are available at most other beamlines at the ALS.



The optical design includes two identical, but oppositely deflecting, toroidal mirrors positioned symmetrically before and after the monochromator. This novel approach simultaneously yields two benefits: high resolution, by providing parallel x-rays for diffraction by the double-crystal monochromator, and a small beam spot, by means of 1:1 focusing of the storage-ring source with a minimum of aberrations. By accepting up to 8 horizontal mrad from an ALS

bend magnet, the beamline delivers >10¹¹ photons/s over most of its photon energy range.

Several interchangeable experimental stations for gas-phase measurements are available for use at the beamline. Two of these instruments use ALS two-bunch operation to do either time-of-flight (TOF) ion mass spectroscopy or photoelectron spectroscopy with five electron-TOF analyzers operating simultaneously. Both devices can use specialized electronics to perform ion-ion or electron-electron coincidence measurements. A third instrument will measure x-ray emission spectroscopy in gases, with both angular and polarization sensitivity. Finally, a small absorption cell is located in the beamline that can be used to perform absorption spectroscopy of gases or thin foils. ■



Approximate photon flux from Beamline 9.3.1. The measurements were taken with Si(111) crystals in the monochromator and at a storage ring current of 200 mA while accepting nearly 3 horizontal mrad from the bend magnet. The decrease observed at high energy is due to the cut-off energy of the beamline mirrors at their nominal angles of incidence. These angles can be varied in situ to shift the cut-off to higher or lower energies.



Absorption spectra. Total yield spectrum of CHCl_3 in the vicinity of the chlorine K edge (left). Total ion yield spectrum of SO_2 in the vicinity of the sulfur K edge (right). All measurements were performed at Beamline 9.3.1 during January/April 2000. Data courtesy of Wayne Stolte (LBNL).

To obtain a proposal form, go to www-als.lbl.gov/als/quickguide/independinvest.html.

For Beamline Information

Fred Schlachter Advanced Light Source Berkeley Lab, MS 7R0222 Berkeley, CA 94720 Tel: (510) 486-4892 Fax: (510) 486-6499 Email: fsschlachter@lbl.gov	Wayne Stolte Advanced Light Source Berkeley Lab, MS 4R0230 Berkeley, CA 94720 Tel: (510) 486-5804 Fax: (510) 495-2111 Email: wcstolte@lbl.gov
--	---

Spokesperson

Dennis W. Lindle
Department of Chemistry
University of Nevada, Las Vegas
Las Vegas, NV 89154-4003
Tel: (702) 895-4426
Fax: (702) 895-4072
Email: lindle@nevada.edu