


# X-Ray Photoemission Electron Microscope (X-PEEM) • Branchline 7.3.1.1

Berkeley Lab • University of California

## Branchline Specifications

Photon Energy Range (eV)	Photon Flux (photons/s/0.1%BW)	Spectral Resolution (E/ΔE)	Polarization	Spatial Resolution (nm)	Samples	Availability
175–1500	$3 \times 10^{12}$ (at 800 eV in 30 × 30 μm spot)	1800 (at 800 eV)	Linear or Circular	20 (NEXAFS) 50 (XMCD)	UHV-Compatible Solids (up to 1 cm <sup>2</sup> in area)	NOW



Schematic layout of Branchline 7.3.1.1.

Beamline 7.3.1 has two branches co-developed with industrial partners. The x-ray photoemission electron microscope (X-PEEM) on Branchline 7.3.1.1 was developed with the IBM Almaden Research Center. It is designed specifically both for near-edge x-ray absorption fine-structure spectroscopy (NEXAFS) and x-ray magnetic circular dichroism (XMCD) and x-ray magnetic linear dichroism (XMLD) spectroscopy of microscopic areas on sample surfaces and for imaging. A separate data sheet describes Branchline 7.3.1.2 for spatially resolved x-ray photoelectron spectroscopy

(microXPS), which operates in parallel with the X-PEEM branch.

The instrument operates over the photon energy range from 175 to 1500 eV with a spectral resolving power between 1000 and 2000 that is adequate for most imaging investigations. A spherical-grating monochromator (SGM) with no entrance slit and a low line-density grating (200 lines/mm) is used to obtain high throughput while still achieving the desired spectral resolution. The radiation source—the center magnet in the triple-bend achromat lattice of the ALS storage ring—is focused to a spot


less than  $30 \times 30 \mu\text{m}$  at the sample. Circularly polarized radiation is obtainable by means of a movable aperture that selects light from above or below the horizontal plane.

The electrostatic optical system of the microscope images secondary electrons at high magnification onto a phosphor and CCD camera by means of a three-element high-voltage (25 kV) objective lens, a transfer lens, an intermediate lens, and a projector lens. The combination of high-voltage operation, a small angle-defining aperture (12  $\mu\text{m}$ ), and stable mechanical construction leads to a resolution of about 20 nm for elemental contrast imaging and 50 nm for XMCD imaging. With electron detection, the samples must be conductive, reasonably flat, UHV-compatible solids.

A flux in the image area greater than  $10^{12}$  photons/s permits real-time studies of elemental, chemical, magnetic, and topographical properties of materials. The yield of electrons over the

entire illuminated area is imaged at high magnification, and, by recording frames sequentially with an incrementally increasing photon energy, one can produce a NEXAFS spectrum for each image point. The same procedure using circularly polarized light yields magnetic images and XMCD spectra and using linearly polarized light yields XMLD images and spectra.

The X-PEEM is equipped with an automated sample-transfer system with a three-sample parking stage and a sample preparation chamber. The sample preparation chamber contains four evaporators, an Ar sputter gun, a sample heater (up to 500 °C), a magnet (up to 1000 Oe), and a LEED system for sample analysis. The endstation is also equipped with an XMCD chamber that allows the acquisition of XMCD spectra in alternating magnetic fields (up to 750 Oe) without spatial resolution (spot size about 1 mm) for reference. ■



*PEEM electron optics layout.*

**This branchline is available to independent investigators with the concurrence of the PRT.**

#### For Branchline Information

Andreas Scholl  
Advanced Light Source  
Berkeley Lab, MS 4-230  
Berkeley, CA 94720  
Tel: (510) 486-4867  
Tel: (510) 495-2111  
Email: a\_scholl@lbl.gov

#### PRT Spokesperson

Joachim Stöhr  
Stanford Synchrotron Radiation  
Laboratory  
MS 69, P.O. Box 4349  
Stanford, CA 94309-0210  
Tel: (650) 926-2570  
Fax: (650) 926-4100  
Email: stohr@ssrl.slac.stanford.edu

#### To Obtain a Proposal Form

Bernie Dixon  
User Services Office Manager  
Advanced Light Source  
Berkeley Lab, MS 6-2100  
Berkeley, CA 94720  
Tel: (510) 486-6722  
Fax: (510) 486-4773  
Email: alsuser@lbl.gov